

## Product Datasheet

### Respiratory Syncytial Virus (RSV) Antibody HM439

<b>Product Name</b>	Respiratory Syncytial Virus (RSV) Antibody
<b>Catalog Number</b>	HM439
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF
<b>Immunogen</b>	Purified native RSV
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Influenza A, Influenza B, and Adenoviruses
<b>Specificity</b>	Fusion protein of RSV strains A & B
<b>Matched Pairs</b>	HM437 (capture)
<b>Preservatives</b>	0.05% Sodium Azide