

Product Datasheet

Troponin I (41-60) Antibody HM456

| | |
|----------------------------|---|
| Product Name | Troponin I (41-60) Antibody |
| Catalog Number | HM456 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, TIA |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified and Filtered Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG1, kappa |
| Clonality | Monoclonal |
| Purity | ≥95% |
| Buffer | 20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl |
| Specificity | Troponin I (epitope 41-60 a.a.) |
| Matched Pairs | HM233 (capture) |