

Product Datasheet

Japanese Encephalitis Virus (JEV-NS1) Antibody HM468

Product Name	Japanese Encephalitis Virus (JEV-NS1) Antibody
Catalog Number	HM468
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Japanese Encephalitis Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with other Flaviviruses
Specificity	Japanese Encephalitis Virus NS1 protein
Preservatives	0.05% Sodium Azide