

## Product Datasheet

### Japanese Encephalitis Virus (JEV-NS1) Antibody HM474

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Japanese Encephalitis Virus (JEV-NS1) Antibody                            |
| <b>Catalog Number</b>      | HM474   |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | LF, ELISA   |
| <b>Immunogen</b>           | Purified Japanese Encephalitis Virus                                      |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Ion Exchange Purified Monoclonal Antibody                                 |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1, kappa   |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | >90% by SDS-PAGE  |
| <b>Buffer</b>              | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride             |
| <b>Cross Reactivity</b>    | No reactivity with other Flaviviruses                                     |
| <b>Specificity</b>         | Japanese Encephalitis Virus NS1 protein                                   |
| <b>Preservatives</b>       | 0.05% Sodium Azide  |