

## Product Datasheet

### Tick-borne Encephalitis Virus (TBEV-NS1) Antibody HM478

<b>Product Name</b>	Tick-borne Encephalitis Virus (TBEV-NS1) Antibody
<b>Catalog Number</b>	HM478
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Purified Tick-borne Encephalitis Virus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	No reactivity with other Flaviviruses
<b>Specificity</b>	Tick-borne Encephalitis Virus NS1 protein
<b>Preservatives</b>	0.05% Sodium Azide