

Product Datasheet

Yellow Fever Virus (YFV-NS1) Antibody HM490

Product Name	Yellow Fever Virus (YFV-NS1) Antibody
Catalog Number	HM490
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Immunogen	Purified Yellow Fever Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with other Flaviviruses
Specificity	Yellow Fever Virus NS1 protein
Preservatives	0.05% Sodium Azide