

Product Datasheet

Anti-h PG II 8102 SPRN-1 100098

Product Name	Anti-h PG II 8102 SPRN-1
Catalog Number	100098
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, IT, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	1.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	KA= 1 x 10 ⁹ 1/M
Associated Products	Recombinant Pepsinogen II antigen, Medix Biochemica 610001 and native Pepsinogen II antigen Lee Biosolutions 440-52
Buffer	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.2-8.0
Cross Reactivity	Human Pepsinogen I 0.9 %, Porcine Pepsinogen < 0.22 %, Gastrin I < 0.22 %, Somatostatin < 0.22 %, Peptide YY < 0.22 %, Bombesin < 0.22 %
Specificity	Antibody recognizes human pepsinogen II
Shelf Life	24 months