

## Product Datasheet

### Methicillin-resistant Staphylococcus aureus (MRSA) Antibody HM504

<b>Product Name</b>	Methicillin-resistant Staphylococcus aureus (MRSA) Antibody
<b>Catalog Number</b>	HM504
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IP, IHC
<b>Immunogen</b>	Purified MRSA protein (mecA)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Monoclonal Antibody
<b>Storage</b>	Lyophilized: 2-8°C (desiccated) Reconstituted: Short Term 2-8°C Long Term: -20°C. Avoid repeated freezing and thawing
<b>Note</b>	Reconstitute with Distilled Water. Recommended starting dilution is 1:5,000.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal