

Product Datasheet

Malaria (pf-HRP-2) Antibody HM506

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM506
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Plasmodium falciparum HRP-2
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	LA278
Preservatives	0.05% Sodium Azide