

Product Datasheet

Pepsinogen I (PGI) Antibody HM510

Product Name	Pepsinogen I (PGI) Antibody
Catalog Number	HM510
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human pepsinogen I
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	PGI
Matched Pairs	HM511 (detection)
Preservatives	0.09% Sodium Azide