

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM513

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM513
Species/Host	Mouse
Tested Applications	LF
Immunogen	Recombinant Dengue Antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Specificity	NS1 protein of Dengue Virus types 1,2,3,4
Matched Pairs	HM932 (capture) & HM931 (detection)
Preservatives	0.1% Proclin 300