

Product Datasheet

Calicivirus, Feline Antibody HM515

Product Name	Calicivirus, Feline Antibody
Catalog Number	HM515
Species/Host	Mouse
Tested Applications	ELISA, WB, IHC, Neutralization, IFA
Immunogen	Feline Calicivirus capsid
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	no reactivity with Feline Herpesvirus, FIV or FeLV
Specificity	Feline Calicivirus VP1 capsid protein
Matched Pairs	HM515 (self-pairs), HM516 (capture) or HM517 (detection)
Preservatives	0.05% Sodium Azide