

## Product Datasheet

### Human Perforin ELISA Kit 850.860

<b>Product Name</b>	Human Perforin ELISA Kit
<b>Catalog Number</b>	850.860
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	62.5-2000 pg/ml
<b>Sensitivity</b>	40 pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with Trail, TRAIL R3, TRAIL R4, Granzyme B, TNF $\alpha$ , IFN $\gamma$ , Fas, Fas L, IL-1b, IL-2
<b>Specificity</b>	Recognizes natural human Perforin
<b>Sample Types</b>	Supernatant, Buffered Solutions, Not suitable for use with serum and plasma.