

Product Datasheet

Parvovirus, Canine Antibody HM543

Product Name	Parvovirus, Canine Antibody
Catalog Number	HM543
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Canine parvovirus
Matched Pairs	HM544 (detection)
Preservatives	0.09% Sodium Azide