

Product Datasheet

Distemper Virus, Canine Antibody HM547

Product Name	Distemper Virus, Canine Antibody
Catalog Number	HM547
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Canine Distemper Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Canine Distemper Virus
Matched Pairs	HM545 or HM546 (capture)
Preservatives	0.09% Sodium Azide