

Product Datasheet

Serum Amyloid A (SAA), Feline Antibody HM558

Product Name	Serum Amyloid A (SAA), Feline Antibody
Catalog Number	HM558
Species/Host	Mouse
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Feline SAA
Matched Pairs	HM558 (self-pairs) or HM559
Preservatives	0.1% Sodium Azide