

## Product Datasheet

### Parainfluenza Virus (CPIV), Canine Antibody HM562

<b>Product Name</b>	Parainfluenza Virus (CPIV), Canine Antibody
<b>Catalog Number</b>	HM562
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IFA
<b>Immunogen</b>	Canine Parainfluenza Virus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG2a, kappa
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	No reactivity with Canine Adenovirus type 2, Canine Coronavirus, Canine Distemper Virus, Canine Parvovirus, FeLV and FIV
<b>Specificity</b>	Canine Parainfluenza Virus
<b>Preservatives</b>	0.05% Sodium Azide