

Product Datasheet

Procalcitonin (PCT) Antibody HM564

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM564
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Procalcitonin
Matched Pairs	HM563 (capture)
Preservatives	0.05% ProClin 300