

## Product Datasheet

### Procalcitonin (PCT) Antibody HM564

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Procalcitonin (PCT) Antibody  |
| <b>Catalog Number</b>      | HM564   |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | LF, ELISA, CLIA   |
| <b>Immunogen</b>           | Recombinant Procalcitonin   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody  |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | >95%  |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.0                                   |
| <b>Specificity</b>         | Procalcitonin   |
| <b>Matched Pairs</b>       | HM563 (capture)   |
| <b>Preservatives</b>       | 0.02% Sodium Azide  |