

Product Datasheet

C-Reactive Protein (CRP) Antibody HM573

Product Name	C-Reactive Protein (CRP) Antibody
Catalog Number	HM573
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Purified human CRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human CRP
Matched Pairs	HM574 (detection)
Preservatives	0.02% Sodium Azide