

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM576

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM576
Species/Host	Mouse
Tested Applications	ELISA, CLIA, CMIA
Immunogen	CK-MB recombinant protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (\leq 1 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	20 mM Phosphate Buffered Saline, pH 7.4
Specificity	Creatine Kinase MB
Molecular Weight	150 kDa
Matched Pairs	HM575 (capture)
Preservatives	None