

Product Datasheet

Coronavirus, Canine Antibody HM591

Product Name	Coronavirus, Canine Antibody
Catalog Number	HM591
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified Canine Coronavirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Canine Coronavirus
Matched Pairs	HM617 (capture) or HM591 (self-pairs)
Preservatives	0.02% Sodium Azide