

Product Datasheet

Coronavirus, Feline Antibody HM592

Product Name	Coronavirus, Feline Antibody
Catalog Number	HM592
Species/Host	Mouse
Tested Applications	ELISA, WB, IHC, IFA
Immunogen	Coronavirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Minimal cross reactivity with Bovine Coronavirus; No reactivity with FeLV, FIV, Feline Calicivirus, Feline Herpes Virus, Canine Adenovirus, Canine Distemper Virus, Canine Parvovirus, Canine Parainfluenza Virus
Specificity	Feline Coronavirus, Canine Coronavirus, Transmissible Gastroenteritis Coronavirus, Ferret Coronavirus
Matched Pairs	HM593 (capture)
Preservatives	0.05% Sodium Azide