

Product Datasheet

Coronavirus, Feline Antibody HM593

Product Name	Coronavirus, Feline Antibody
Catalog Number	HM593
Species/Host	Mouse
Tested Applications	ELISA, WB, IHC, IFA
Immunogen	Coronavirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Canine Adenovirus, Canine Distemper Virus, Canine Parvovirus, Canine Parainfluenza Virus
Specificity	Feline Coronavirus, Canine Coronavirus, Transmissible Gastroenteritis Coronavirus
Matched Pairs	HM592 (detection)
Preservatives	0.05% Sodium Azide