

Product Datasheet

UR-144 Metabolite Antibody HM596

Product Name	UR-144 Metabolite Antibody
Catalog Number	HM596
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	UR-144-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	UR-144 5%
Specificity	UR-144 5-hydroxypentyl, UR-144 5-pentanoic acid, UR-144 4-hydroxypentyl
Matched Pairs	LA385 (capture)
Preservatives	0.05% Sodium Azide