

Product Datasheet

Insulin Antibody HM604

Product Name	Insulin Antibody
Catalog Number	HM604
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Recombinant Insulin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-70°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	≥90%
Buffer	15 mM Phosphate Buffer, pH 7.2 with 0.15M NaCl
Matched Pairs	HM603 (capture)
Preservatives	0.09% Sodium Azide