

## Product Datasheet

### Insulin Antibody HM604

<b>Product Name</b>	Insulin Antibody
<b>Catalog Number</b>	HM604
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Insulin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	-70°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥90%
<b>Buffer</b>	15 mM Phosphate Buffer, pH 7.2 with 0.15M NaCl
<b>Matched Pairs</b>	HM603 (capture)
<b>Preservatives</b>	0.09% Sodium Azide