

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM606

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM606
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	AMH
Matched Pairs	HM605 (capture)
Preservatives	0.05% Proclin 300