

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM607

Product Name Anti-Mullerian Hormone (AMH) Antibody

Catalog Number HM607

Species/Host Mouse

Tested Applications ELISA, CLIA

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage 2-8°C.

Clonality Monoclonal

Purity ≥90%

Buffer Phosphate Buffered Saline, pH 7.4

Matched Pairs HM608 (detection)

Preservatives 0.05% ProClin 300