

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM608

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM608
Species/Host	Mouse
Tested Applications	ELISA, CLIA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Limit of Detection for CLIA \leq 0.05 ng/ml, IFA 0.03 ng/ml
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	AMH N-terminal domain
Matched Pairs	HM607 (capture)
Preservatives	None