

Product Datasheet

Pinaca Antibody HM614

| | |
|----------------------------|---|
| Product Name | Pinaca Antibody |
| Catalog Number | HM614 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | AB-PINACA, AB-PINACA 4-hydroxypentyl & AB-PINACA 5-hydroxypentyl 100%; AB-PINACA 5-pentanoic acid 95%; APINACA & APINACA 5-hydroxypentyl 1% |
| Specificity | AB-PINACA |
| Matched Pairs | LA366 (capture) |
| Preservatives | 0.05% Sodium Azide |