

Product Datasheet

Pinaca Antibody HM614

Product Name	Pinaca Antibody
Catalog Number	HM614
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	AB-PINACA, AB-PINACA 4-hydroxypentyl & AB-PINACA 5-hydroxypentyl 100%; AB-PINACA 5-pentanoic acid 95%; APINACA & APINACA 5-hydroxypentyl 1%
Specificity	AB-PINACA
Matched Pairs	LA366 (capture)
Preservatives	0.05% Sodium Azide