

Product Datasheet

Pinaca Antibody HM615

Product Name Pinaca Antibody

Catalog Number HM615

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Buffer 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity AB-PINACA, AB-PINACA 4-hydroxypentyl & AB-PINACA 5-hydroxypentyl 100%; AB-PINACA 5-pentanoic acid 95%; APINACA & APINACA 5-hydroxypentyl 1%

Specificity AB-PINACA

Matched Pairs LA366 (capture)

Preservatives 0.05% Sodium Azide