

## Product Datasheet

### Pinaca Antibody HM616

**Product Name** Pinaca Antibody

**Catalog Number** HM616

**Species/Host** Mouse

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Purified Monoclonal Antibody

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG1

**Clonality** Monoclonal

**Buffer** 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

**Cross Reactivity** APINACA 5-hydroxypentyl 100%; APINACA 80%; UR-144 4-hydroxypentyl 10%; UR-144 5-hydroxypentyl; UR-144, UR-144 5-pentanoic acid, AB-PINACA, AB-PINACA 5-hydroxypentyl & AB-PINACA 5-pentanoic acid  $\leq$  3%

**Specificity** APINACA

**Matched Pairs** LA366 (capture)

**Preservatives** 0.05% Sodium Azide