

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM618

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM618
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Limit of Detection for Lateral Flow 1 ng/ml (colloidal gold) and \leq 20 ng/ml fluorescent label
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	AMH N-terminal domain
Matched Pairs	HM619 (detection)
Preservatives	0.02% Proclin 300