

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM619

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM619
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	AMH N-terminal domain
Matched Pairs	HM618 (capture)
Preservatives	0.05% ProClin 300