

Product Datasheet

Usutu Virus (USUV-NS1) Antibody HM624

Product Name	Usutu Virus (USUV-NS1) Antibody
Catalog Number	HM624
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Antigen (Vienna 2001)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No cross-reactivity with other Flaviviruses
Specificity	Usutu Virus
Matched Pairs	HM622 or HM623 (capture)
Preservatives	0.05% Sodium Azide