

Product Datasheet

Mycobacterium tuberculosis (TB) Antibody HM627

Product Name	Mycobacterium tuberculosis (TB) Antibody
Catalog Number	HM627
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Mycobacterium tuberculosis
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	M. tuberculosis BlaC
Preservatives	0.05% Sodium Azide