

Product Datasheet

Fentanyl Antibody HM628

Product Name	Fentanyl Antibody
Catalog Number	HM628
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Fentanyl-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Acetyl Fentanyl 100%; Norfentanyl & Acetyl Norfentanyl 90%; Alfentanil 30%; Risperidone 20%; Carfentanil 5%
Specificity	Fentanyl
Matched Pairs	LA395 (capture)
Preservatives	0.05% Sodium Azide