

## Product Datasheet

### Rotavirus (VP6) Antibody HM630

<b>Product Name</b>	Rotavirus (VP6) Antibody
<b>Catalog Number</b>	HM630
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Recombinant Rotavirus VP6 (LA397)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Campylobacter, Clostridium difficile, E. coli, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica, Enterovirus, Norovirus, Astrovirus, Adenovirus, Cryosporidium parvum, Giardia lamblia, H. pylori, Entamoeba histolytica, S. aureus
<b>Specificity</b>	Rotavirus VP6
<b>Matched Pairs</b>	HM630 (self-pairs)
<b>Preservatives</b>	0.09% Sodium Azide