

Product Datasheet

Calprotectin Antibody HM632

Product Name	Calprotectin Antibody
Catalog Number	HM632
Species/Host	Mouse
Tested Applications	LF, ELISA, TIA
Immunogen	Recombinant Calprotectin (LA399)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with hemaglobin (human, bovine and porcine), transferrin (human, bovine and porcine) and lactoferrin (human and bovine)
Specificity	Human Calprotectin S100A9 (native and recombinant)
Matched Pairs	HM632 (self-pairs against calprotectin tetramer)
Preservatives	0.09% Sodium Azide