

Product Datasheet

Mycoplasma pneumoniae Antibody HM633

Product Name	Mycoplasma pneumoniae Antibody
Catalog Number	HM633
Species/Host	Mouse
Tested Applications	LF
Immunogen	M. pneumoniae FH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with M. fermentans, M. genitalium, M. salvarium, M. hominis and M. orale
Specificity	M. pneumoniae P30
Matched Pairs	HM635 (detection) and HM634 or HM636 (capture)
Preservatives	0.05% Sodium Azide