

## Product Datasheet

### Human MICA ELISA Kit 850.910

<b>Product Name</b>	Human MICA ELISA Kit
<b>Catalog Number</b>	850.910
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	156.25-5000 pg/ml
<b>Sensitivity</b>	123 pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with IL-7, IL-12p40, IL-2R, Trail, CD31, IL-13R $\alpha$ 2, IL-5, gp130, TNF $\alpha$ .
<b>Specificity</b>	Recognizes both natural and recombinant human MICA
<b>Sample Types</b>	Supernatant, Buffered Solutions. Not suitable for use with serum and plasma.