

Product Datasheet

Mycoplasma pneumoniae Antibody HM636

| | |
|----------------------------|--|
| Product Name | Mycoplasma pneumoniae Antibody |
| Catalog Number | HM636 |
| Species/Host | Mouse |
| Tested Applications | LF |
| Immunogen | M. pneumoniae FH |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >95% by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | No reactivity with M. fermentans, M. genitalium, M. salvarium, M. hominis and M. orale |
| Specificity | M. pneumoniae P30 |
| Matched Pairs | HM633 (detection) |
| Preservatives | 0.05% Sodium Azide |