

## Product Datasheet

### Adenovirus Type 1 & 2, Canine Antibody HM641

<b>Product Name</b>	Adenovirus Type 1 & 2, Canine Antibody
<b>Catalog Number</b>	HM641
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IHC, IFA
<b>Immunogen</b>	Canine Adenovirus Types 1 & 2
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	No reactivity with Canine Coronavirus, Canine Distemper Virus, Canine Parvovirus, Canine Parainfluenza Virus
<b>Specificity</b>	Canine Adenovirus Types 1 & 2
<b>Matched Pairs</b>	HM642 (capture or detection)
<b>Preservatives</b>	0.05% Sodium Azide