

Product Datasheet

Adenovirus Type 1 & 2, Canine Antibody HM642

Product Name	Adenovirus Type 1 & 2, Canine Antibody
Catalog Number	HM642
Species/Host	Mouse
Tested Applications	ELISA, IFA
Immunogen	Canine Adenovirus Types 1 & 2
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Canine Coronavirus, Canine Distemper Virus, Canine Parvovirus, Canine Parainfluenza Virus
Specificity	Canine Adenovirus Types 1 & 2
Matched Pairs	HM641 (capture or detection)
Preservatives	0.05% Sodium Azide