

Product Datasheet

Aldosterone Antibody HM651

Product Name	Aldosterone Antibody
Catalog Number	HM651
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Aldosterone-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	17 α -Hydroxyprogesterone, Estrone, 17 α -Ethinylestradiol, 17 β -Estradiol, Estriol, Testosterone, Testosterone glucuronide, 5 α -Dihydrotestosterone, Epi-testosterone, 6 β -Hydroxytestosterone, 17 α -Methyltestosterone, Progesterone, Androstenedione <1%
Specificity	Aldosterone
Matched Pairs	LA427 or LA468 (capture)
Preservatives	0.05% Sodium Azide