

Product Datasheet

Aldosterone Antibody HM653

Product Name Aldosterone Antibody

Catalog Number HM653

Species/Host Mouse

Tested Applications LF, ELISA

Immunogen Aldosterone-KLH

Brand EastCoast Bio

Form/Appearance Ion Exchange Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2b

Clonality Monoclonal

Purity >90% by SDS-PAGE

Buffer 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity 17 α -Hydroxyprogesterone, Estrone, 17 α -Ethinylestradiol, 17 β -Estradiol, Estriol, Testosterone, Testosterone glucuronide, 5 α -Dihydrotestosterone, Epitestosterone, 6 β -Hydroxytestosterone, 17 α -Methyltestosterone, Progesterone, Androstenedione <1%

Specificity Aldosterone

Matched Pairs LA427 or LA468 (capture)

Preservatives 0.05% Sodium Azide