

Product Datasheet

Estrone-3-Glucuronide (E1G) Antibody HM658

Product Name	Estrone-3-Glucuronide (E1G) Antibody
Catalog Number	HM658
Species/Host	Rat
Tested Applications	LF, ELISA
Immunogen	Estrone-3-Glucuronide-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Estriol <6%; Estriol-3-glucuronide, Estradiol-3-glucuronide <5%; 17 α -Ethinylestradiol <3%; 17 α -hydroxyprogesterone, Testosterone, Progesterone, Pregnenediol-3-glucuronide, Aldosterone, Androstenedione, Estrone sulfate <1%
Specificity	Estrone and Estrone-3-glucuronide
Matched Pairs	LA432 (capture)
Preservatives	0.05% Sodium Azide