

## Product Datasheet

### Serum Amyloid A (SAA), Feline Antibody HM702

<b>Product Name</b>	Serum Amyloid A (SAA), Feline Antibody
<b>Catalog Number</b>	HM702
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Feline Serum Amyloid A
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	sensitivity $\leq$ 10 ng/ml
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	10 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride
<b>Specificity</b>	Feline Serum Amyloid A
<b>Matched Pairs</b>	HM703 (capture)
<b>Preservatives</b>	0.1% Sodium Azide